

## CELL AND TISSUE CULTURE DISHES CELLATTACH® SURFACE TREATED



**37.58€ - 46.68€**

Cell and Tissue Culture Dishes are ideal containers for cell growth and yields aim on little and medium input volume. Surface Treated CellATTACH. 35 mm of diameter, 8,5 cm<sup>2</sup> of cell grown area and Sterile, or 60 mm of diameter, 21,2 cm<sup>2</sup> of cell grown area and Sterile. Box of 10 bags x 24 units.

**SKU:** CAD010060 | **Category:** [Cell and Tissue Culture Dishes](#) |

### VARIATIONS

Image	SKU	Price	Description	Unidades
	CAD010035	37.58€		35 mm (204 units)
	CAD010060	46.68€		60 mm (240 units)

PRODUCT DESCRIPTION

Cell and Tissue Culture Dishes are ideal containers for cell growth and yields aim on little and medium input volume. Surface Treated CellATTACH. 35 mm of diameter, 8,5 cm2 of cell grown area and Sterile, or 60 mm of diameter, 21,2 cm2 of cell grown area and Sterile

Flat bottom uniform wall thickness ensures distortion-free bottom.

Dish surface is smooth and free from striation to maximize usable area for growth.

The rim on upper side of the lid mates with the dish brim for easy and secure stack.

Lids with several little chimbs to shape vents are available for very effective gas exchange.

Sterilized by gamma irradiation.

Non-pyrogenic.

ADDITIONAL INFORMATION

Weight	4.5 kg
Dimensions	40.7 × 27.6 × 21.5 cm
Unidades	35 mm (204 units), 60 mm (240 units)



