


## CELL AND TISSUE CULTURE FLASK STANDARD SURFACE TREATED AND STANDARD CAP



**51.79€ - 54.20€**


Cell and Tissue Culture Flasks are perfect for cell growth and cell yields aim on little and medium input volume. Surface treated standard and standard cap, sterile. Volumes: 25 Appro. Cell growth area 12,5 cm<sup>2</sup>, or 50 ml Appro. Cell growth area 25 cm<sup>2</sup>. Packaging: 200 units

**SKU:** TCF011050 | **Category:** [Cell and Tissue Culture Flask](#) |

### VARIATIONS

Image	SKU	Price	Description	Unidades
	TCF011025	51.79€		25 ml (200 units)
	TCF011050	54.20€		50 ml (200 units)

GALLERY IMAGES



PRODUCT DESCRIPTION

Cell and Tissue Culture Flasks are perfect for cell growth and cell yields aim on little and medium input volume. Surface treated standard and standard cap, sterile. Volumes: 25 Appro. Cell growth area 12,5 cm2, or 50 ml Appro. Cell growth area 25 cm2.  
Flat bottom uniform wall thickness ensures distortion-free bottom.

Dish surface is smooth and free from striation to maximize usable area for growth.

The rim on upper side of the lid mates with the dish brim for easy and secure stack.

Lids with several little chimbs to shape vents are available for very effective gas exchange.

Sterilized by gamma irradiation.

Non-pyrogenic.

**ADDITIONAL INFORMATION**

<b>Weight</b>	5.40 kg
<b>Dimensions</b>	52.0 × 32.0 × 24.3 cm
<b>Unidades</b>	25 ml (200 units), 50 ml (200 units)



