


CELL AND TISSUE CULTURE PLATES STANDARD, SURFACE TREATED, U BOTTOM, 96 WELLS



74.68€

Cell and Tissue culture plates Standard, Surface treated, U bottom, 96 wells. Max volume single well 0.33 ml. Sterilized and non-pyrogenic. Packaging: 100 units

SKU: TCP012096 | **Category:** [Cell and Tissue Culture Plates](#) |

GALLERY IMAGES



PRODUCT DESCRIPTION

Cell and Tissue culture plates Standard, Surface treated, U bottom, 96 wells. Max volume single well 0.33 ml. Sterilized and non-pyrogenic.

Cell and Tissue Culture Plates are ideal for cell growth and cell yields on multiple, compare and other analyse.

Uniform well volume ensures equal growth surface area

Flat well bottom and round bottom plate.

Well surface is smooth and free from striation to maximize usable growth area.

Raised rims on wells with the uniform rings on the lid to reduce evaporation.

Single position lid reduces the risks of cross-contamination and the handling mistakes.

Wells are labeled with alphanumeric code for easy identification.

Suitable for all common instruments and automation.

Sterilized by gamma irradiation.

1

