


CELL AND TISSUE CULTURE PLATES STANDARD, SURFACE TREATED








60.21€ - 67.45€

Cell and Tissue culture plates Standard, Surface treated, Available: 6, 12, 24, 48 or 96 (flat bottom) wells. Max volume single well 17.0 ml, 6.80 ml, 3.50 ml, 1.55 ml or 0.39 ml respectively. Sterilized and non-pyrogenic. Packaging: 100 units

SKU: TCP011006 | **Category:** [Cell and Tissue Culture Plates](#) |

VARIATIONS

Image	SKU	Price	Description	Unidades
	TCP011012	62.63€		12 wells (100 units)
	TCP011024	62.63€		24 wells (100 units)
	TCP011048	67.45€		48 wells (100 units)
	TCP011006	60.21€		6 wells (100 units)
	TCP011096	60.21€		96 wells, flat bottom (100 units)

GALLERY IMAGES



PRODUCT DESCRIPTION

Cell and Tissue culture plates Standard, Surface treated. Available: 6, 12, 24, 48 or 96 (flat bottom) wells. Max volume single well 17.0 ml, 6.80 ml, 3.50 ml, 1.55 ml or 0.39 ml respectively. Sterilized and non-pyrogenic.

Cell and Tissue Culture Plates are ideal for cell growth and cell yields on multiple, compare and other analyse.

Uniform well volume ensures equal growth surface area

Flat well bottom and round bottom plate.

Well surface is smooth and free from striation to maximize usable growth area.

Raised rims on wells with the uniform rings on the lid to reduce evaporation.

Single position lid reduces the risks of cross-contamination and the handling mistakes.

Wells are labeled with alphanumeric code for easy identification.

Suitable for all common instruments and automation.

Sterilized by gamma irradiation.

ADDITIONAL INFORMATION

Weight	6.55 kg
Dimensions	59.2 × 40.2 × 31.4 cm
Unidades	12 wells (100 units), 24 wells (100 units), 48 wells (100 units), 6 wells (100 units), 96 wells, flat bottom (100 units)

