

## CELLBOX SURFACE TREATED AND STERILE



**69.07€ - 174.57€**


CellBOX 1 layer or up to 10 layers, Surface Treated and Sterile. Dnasa & Rnasa free and Non-pyrogenic. Packaging: 6 units

**SKU:** CBX010001 | **Category:** [CellBOX™](#) |

## VARIATIONS

Image	SKU	Price	Description	Unidades
	CBX010001	69.07€		1 layer (10 units)
	CBX110010	174.57€		up to 10 layers (10 units)

GALLERY IMAGES



PRODUCT DESCRIPTION

The culture device with monolayer or up to 10 layers, can realize large-scale cell culture by one time.  
Made from high grade gamma resistance polystyrene.  
The culture device with monolayer or up to 10 layers can realize large-scale cell culture by one time.  
Easy to use and standardized operation ensures uniform cell growth.  
Suitable for most adherent cell lines.  
Linear amplification, easy to expand to the industrial scale.

Transparent material can be observed directly under a microscope.

Modular design can be multi-level series or parallel.

It can be equipped with a variety of online monitoring instruments.

Strict integrity tested

Sterilized by gamma irradiation

DNase/RNase-free & Non-pyrogenic

ADDITIONAL INFORMATION

Weight	0.32 kg
Dimensions	30.7 × 16.7 × 22.9 cm
Unidades	1 layer (10 units), up to 10 layers (10 units)



