

## **CENTRIFUGE TUBE, CONICAL BOTTOM AND EASYFLIP TOP CAP, RE-SEALABLE BAG, STERILIZED**



**44.10€ - 60.54€**

Centrifuge tube conical bottom and EasyFlip top cap. Volumes: 15 or 50 ml. Sterilized.  
Packaging: Box with 20 re-sealable bags x 25 units.

**SKU:** CFT211150 | **Category:** [Centrifuge Tubes](#) |

### **VARIATIONS**

Image	SKU	Price	Description	Unidades
	CFT211150	44.10€		15 ml (re-sealable bag)
	CFT221150	49.53€		15 ml (Paper rack)
	CFT211500	60.54€		50 ml (re-sealable bag)

**PRODUCT DESCRIPTION**

Centrifuge tube conical bottom, Flip top cap and sterilized. Volumes: 15 or 50 ml.  
Made from high grade gamma resistance polypropylene  
One hand easy flip to open the cap  
Max.RCF: 9,400xg  
Black printed graduations and white writing area can withstand chloroform  
Autoclavable at 121°C and freezable to -80°C  
Packaged in zip sealed bags or environmental friendly, recyclable paper racks  
Leak proof  
Gamma irradiation sterilized

**ADDITIONAL INFORMATION**

<b>Weight</b>	4.74 kg
<b>Dimensions</b>	0.52 × 0.32 × 0.28 cm
<b>Unidades</b>	15 ml (Paper rack), 15 ml (re-sealable bag), 50 ml (re-sealable bag)



