

CENTRIFUGE TUBE, CONICAL BOTTOM AND FLAT CAP, 15 ML, STERILIZED





34.73€ - 36.61€

Centrifuge tube conical bottom and flat cap. Volume 15 ml. Sterilized, DNase/RNase free and Non-pyrogenic. Re-sellable bag or in rack. Packaging: Box with 500 units (20 bags x 25 units)

SKU: CFT011150 | **Category:** [Centrifuge Tubes](#) |

VARIATIONS

Image	SKU	Price	Description	Unidades
	CFT011150	36.61€		15 ml (Bulk)
	CFT021150	34.73€		15 ml (Rack)

GALLERY IMAGES



PRODUCT DESCRIPTION

Centrifuge tube conical bottom and flat cap. Volume 15 ml. Sterilized, DNase/RNase free and Non-pyrogenic. presentation: in rack or bulk.

Longer length screw caps with sealing ring prevent any leakage.

Easy-to-read black graduations are accurate to $\pm 2\%$, 1ml increments (15ml) or 2.5ml increments (50ml).

With a large unerasable frosting white printed writing area.

Both the graduations and writing areas are chloroform-resistant.

Engraved graduation at the conical bottom of each tube.

Max rotate speed up to 12,000xg for conical bottom tubes, while 6,000xg for self-standing tubes.

Autoclavable at 121°C and freezable to -80°C.

Leak-proof

ADDITIONAL INFORMATION

Weight	4.22 kg
Dimensions	51.5 × 31.3 × 19.2 cm
Unidades	15 ml (Bulk), 15 ml (Rack)

