

CENTRIFUGE TUBE, SELF-STANDING BOTTOM AND FLAT CAP, 50 ML, STERILIZED



52.97€

Centrifuge tube self-standing bottom and flat cap. Volume 50 ml. Sterilized, DNase/RNase free and Non-pyrogenic. Packaging: Box with 20 Re-sealable bags x 25 units.

SKU: CFT111500 | **Category:** [Centrifuge Tubes](#) |

PRODUCT DESCRIPTION

Centrifuge tube self-standing bottom and flat cap. Volume 50 ml. Sterilized, DNase/RNase free and Non-pyrogenic.

Longer length screw caps with sealing ring prevent any leakage.

Easy-to-read black graduations are accurate to $\pm 2\%$, 1ml increments (15ml) or 2.5ml increments (50ml).

With a large unerasable frosting white printed writing area.

Both the graduations and writing areas are chloroform-resistant.

Engraved graduation at the conical bottom of each tube.

Max rotate speed up to 12,000xg for conical bottom tubes, while 6,000xg for self-standing tubes.

Autoclavable at 121°C and freezable to -80°C.

Leak-proof

ADDITIONAL INFORMATION

Weight	8.32 kg
Dimensions	57.8 × 52.7 × 28.5 cm
Unidades	

