

DISPOSABLE GREEN L SHAPED CELL SPREADER. PP. STERILE





7.25€ - 15.70€

Disposable green L-shaped cell spreaders are ideal for spreading cell or bacterial cultures evenly. Available: 100 units in individually Wrapped or 500 units packaged 10 units per pack.

SKU: CSP012014 | **Category:** [Cell Spreader](#) |

VARIATIONS

Image	SKU	Price	Description	Unidades
	CSP011014	7.25€		Individually Wrapped (100 units)
	CSP012014	15.70€		10 units per pack (500 units)

GALLERY IMAGES



PRODUCT DESCRIPTION

Disposable green L-shaped cell spreaders are ideal for spreading cell or bacterial cultures evenly. The L-shaped allows application on the entire dish or plate surface. Smooth spreading surface with slight upward return reduces the risk of damaging agar sample. Spreaders are manufactured of polypropylene and pre-sterilized, eliminating the need for flaming or autoclaving. Sterilized by gamma irradiation

Non-pyrogenic

ADDITIONAL INFORMATION

Weight	1.40 kg
Dimensions	26.0 × 21.6 × 10.0 cm
Unidades	10 units per pack (500 units), Individually Wrapped (100 units)

