


## ELISA PLATE, FIXED FLAT BOTTOM, 96 WELLS




**76.68€**

ELISA Plate, fixed flat bottom, 96 wells. Medium Binding capacity and sterilized. Packaging: 200 units

**SKU:** FEP200096 | **Category:** [Elisa Plates](#) |

## GALLERY IMAGES



## PRODUCT DESCRIPTION

ELISA Plates are optimal products for ELISA and provide reliable performance in binding assays when consistent coating of wells is required.

Well surface is uniform, smooth and free from striation to eliminate error.

CV of transmittance is less than 5.00%.

Wells are labeled with alphanumeric code for easy identification.

Sterilized by gamma irradiation.

1



