

FILTER UPPER CUP, CA 0,22 MM, 500 ML, NON-PYROGENIC.



Filter Upper Cup, CA, membrane pore size 0,22 µm, volume 500 ml. Diameter 75 mm. Non-pyrogenic. Packaging: 24 units

SKU: FCA216500 | **Category:** [Filter Upper Cup](#) |

PRODUCT DESCRIPTION

Filter Upper Cup, membrane pore size 0,22 µm, volume 500 ml. Diameter 75 mm.
Membrane CA.
Light weight and heavy wall construction.

Designed wide and easy access bottle mouth for efficiently and stably pour out.
Non-pyrogenic.

1

THERE ARE NO REVIEWS YET.

