

MICRO CENTRIFUGE TUBE, CONICAL, 1,5 ML. EASYFLIP. NON-STERILIZED



5.07€

Micro Centrifuge tube conical bottom and EasyFlip. Volume 1,5 ml. Non-sterilized, DNase/RNase free and Non-pyrogenic. Packaging: 500 units

SKU: CFT000015 | **Category:** [Micro Centrifuge Tubes](#) |

PRODUCT DESCRIPTION

Micro Centrifuge tube conical bottom and EasyFlip. Volume 1,5 ml. Non-sterilized.
Made from high clarity and high grade gamma resistance polypropylene.
One hand easy flip to open the cap.

Engraved graduation ensure accuracy.
Frosted body surface provide easy and legible mark.
DNase/RNase free and non-pyrogenic.
Suitable for all standard racks and rotors.
Boil proof.
Maximum rotate speed up to 25,000xg.
Autoclavable at 121°C and freezable to -80°C.

ADDITIONAL INFORMATION

| | |
|------------|-----------------------|
| Weight | 0.56 kg |
| Dimensions | 6.98 × 5.83 × 2.13 cm |
| Unidades | |

