

## MICRO CENTRIFUGE TUBE, CONICAL AND SCEW CAP. 2 ML, STERILIZED



**22.85€**

Micro Centrifuge tube conical bottom. Volume 1,5 ml. Non-sterilized or sterile, DNase/RNase free and Non-pyrogenic

**SKU:** CFT003020 | **Category:** [Micro Centrifuge Tubes](#) |

### PRODUCT DESCRIPTION

Micro Centrifuge tube conical bottom and scREW cap. Volume 2 ml. Sterilized.  
Engraved graduation ensure accuracy.  
Flat and frost caps surface together with smooth and frosted body surface.

DNase/RNase free and non-pyrogenic.  
Autoclavable at 121°C and freezable to -80°C

1

**THERE ARE NO REVIEWS YET.**

