

MICRO CENTRIFUGE TUBE, CONICAL BOTTOM, 1,5 ML, EASYFLIP





6.09€ - 6.34€

Micro Centrifuge tube conical bottom. Volume 1,5 ml. Non-sterilized or sterile, DNase/RNase free and Non-pyrogenic. Packaging: 500 units

SKU: CFT002015 | **Category:** [Micro Centrifuge Tubes](#) |

VARIATIONS

Image	SKU	Price	Description	Unidades
	CFT-003-015	6.34€		Sterile (500 units)
	CFT002015	6.09€		Non-esterilized (500 units)

PRODUCT DESCRIPTION

Micro Centrifuge tube conical bottom. Volume 1,5 ml. Sterile or Non-sterilized.
Engraved graduation ensure accuracy.
Flat and frost caps surface together with smooth anf frosted body surface.
DNase/RNase free and non-pyrogenic.
Autoclavable at 121°C and freezable to -80°C

ADDITIONAL INFORMATION

Weight	0.56 kg
Dimensions	6.98 × 5.83 × 2.13 cm
Unidades	Non-esterilized (500 units), Sterile (500 units)

