

## MICRO CENTRIFUGE TUBE WITH LID, 1,5 ML, NON-STERILIZED



**206.48€**

Micro centrifuge round bottom tube with lid and maximum rotate speed up to 14.000 g. Volume 1,5 ml, Non-sterilized. Packaging: 200 units

**SKU:** CFT108015 | **Category:** [Centrifuge Tubes](#) |

### PRODUCT DESCRIPTION

Micro centrifuge round bottom tube of Polypropylene, with lid and maximum rotate speed up to 14.000 g. Volume 1,5 ml.

Non-sterilized. Re-sealable bag. Dnasa / Rnasa free and Non-pyrogenic

ADDITIONAL INFORMATION

Weight	6.31 kg
Dimensions	0.43 × 0.38 × 0.41 cm
Unidades	



