

MICRO VOLUME PIPETTORs



80.00€

Variable volume single-channel pipette. Volumes: 0,1-2,5 µl, 20-200 µl and 100-1000 µl.
Packaging: Unit

SKU: N/A | **Category:** [Micro Volume Pipettors](#) |

VARIATIONS

Image	SKU	Price	Description	Unidades
	SPA200125	80.00€		0,1-2,5 µl
	SPA211000	80.00€		100-1000 µl
	SPA220200	80.00€		20-200 µl

GALLERY IMAGES





PRODUCT DESCRIPTION

Variable volume single-channel pipette. Volumes: 0,1-2,5 µl, 20-200 µl and 100-1000 µl

Lightweight format and resistant to UV light.

Fully autoclavable.

Simple recalibration mechanism.

Minimum resistance in the piston and the tip ejector.

Includes calibration tool.

It has a color-coded push button to facilitate volume recognition.

ADDITIONAL INFORMATION

Weight	4.80 kg
Dimensions	0.30 × 0.20 × 0.10 cm
Unidades	0,1-2,5 µl, 0.5-10 uL, 10-100 uL, 100-1000 µl, 2-20 uL, 20-200 µl, 5-50 uL

