

## MULTICHANNEL MICRO PIPETTOR, 8 CHANNELS







**109.88€**

Multichannel Micro Pipettor, 8 channels. Volumes: 0.5-10 µL, 5-50 µL, 10-100 µL and 20-200 µL.  
Packaging: unit

**SKU:** SPA008100 | **Category:** [Multichannel Micro Pipettors](#) |

## VARIATIONS

Image	SKU	Price	Description	Unidades
	SPA008010	109.88€		0,5-10µl
	SPA008100	109.88€		10-100µl
	SPA008200	109.88€		20-200µl
	SPA008050	109.88€		5-50µl

## PRODUCT DESCRIPTION

Automatic 8-channel Jetpip pipette. Volumes: 0.5-10  $\mu$ L, 5-50  $\mu$ L, 10-100  $\mu$ L and 20-200  $\mu$ L.

Soft and smart TPE grip.

Easily removable and fully autoclavable collector. With 360° turn option to facilitate operation with the right or left hand.

Soft tip loading and ejection. Ejection of the tip with one hand.

With ring at the end that prevents the leakage of liquids.

Color coded for better identification

Compatible with most types of tips on the market.

Calibration report included in each pipette.

ADDITIONAL INFORMATION

Weight	0.61 kg
Dimensions	0.29 × 0.14 × 0.05 cm
Unidades	0,5-10µl, 10-100µl, 20-200µl, 5-50µl



