

PCR TUBE WITH FLAT CAP, 0,2 ML



26.63€ - 357.54€

PCR tube with flat cap, 0,2 ml, single or 8 strips. Packaging: bag of 1000 units single or 10 bags of 125 strip of 8 tubes

SKU: PCR000200 | **Category:** [PCR tubes](#) |

VARIATIONS

Image	SKU	Price	Description	Unidades
	PCR000200	26.63€		single (1000 units)
	PCR100200	357.54€		8 strips (125 strips)

PRODUCT DESCRIPTION

PCR tube with flat cap, 0,2 ml, single or 8 strips

PCR products are manufactured, from prime virgin, polypropylenes. This result in tubes, strips and, tips that exhibit a perfect balance between transparency, softness, robustness, antistatic characteristics and gas tightness.

Manufactured from prime virgin polypropylene.
Ultra thin wall design for efficient heat transfer.
Compatible with standard 96-well heat blocks.
High light transmittance.
Sterilized.
DNase/RNase-free.

Package in easy opening plastic bag and white box.
Every case has printed lot No. for quality traceability.

ADDITIONAL INFORMATION

Weight	2.50 kg
Dimensions	34.6 × 27.5 × 30.0 cm
Unidades	8 strips (125 strips), single (1000 units)

