

## PIPETTE CONTROLLER FOR LIQUID HANDLING, VOLUME 0.1-100 ML, 0,20 MM



**107.77€**

Automatic pipettor for liquid handling with volumes of 0.1-100 ml, 0,20  $\mu$ m. Universal nozzle compatible with most pipettes on the market. Packaging: Unit

**SKU:** SPA004220 | **Category:** [Pipette Controller](#) |

### PRODUCT DESCRIPTION

Automatic pipettor for liquid handling with volumes of 0.1-100 ml, 0,20  $\mu$ m. Universal nozzle compatible with most pipettes on the market.

Lightweight, ergonomic and safe automatic pipettor, used correctly, prevents contact with suction liquids. It has a universal nozzle, compatible with most pipettes on the market.

Speed adjustment made easy with your fingertips.  
Lightweight, balanced, ergonomic design allows easy pipetting.  
Speed settings.  
Rechargeable lithium polymer wireless battery with long operating time.  
Smooth adjustment of pump speed.  
Compatible with most 0.1-100 ml glass and plastic pipettes.  
Filling capacity of a 25 ml pipette in <5 seconds.  
Quick release of the suction cone for easy exchange of membrane filters.

Included:  
1 charger  
1 wall bracket  
2 membrane filters  
1 desktop stand

ADDITIONAL INFORMATION

Weight	22.5 kg
Dimensions	58.0 × 38.0 × 50.0 cm
Unidades	

**THERE ARE NO REVIEWS YET.**

