

## POLYCARBONATE (PC) MEMBRANE TISSUE CULTURE PLATE INSERTS, 12 WELLS




**69.00€ - 75.90€**

Semitransparent polycarbonate membrane with 12 wells and a pore size of 0.1 or 1.0  $\mu\text{m}$ .  
Sterilized. Packaging: 48 units (box with 4 packs of 12 units)

**SKU:** TCS000012 | **Category:** [Tissue Culture Plate Insert](#) |

### VARIATIONS

Image	SKU	Price	Description	Unidades
	TCS000012	75.90€		pore size 0.1µm (48 units)
	TCS005012	69.00€		pore size 1.0 µm (48 units)

GALLERY IMAGES



PRODUCT DESCRIPTION

The polycarbonate plate inserts feature a thin, semitransparent polycarbonate membrane, with pore size 0.1µm. All inserts are treated for optimal cell attachment. The polycarbonate plate inserts are sterile and assembled with well plates and resist most fixing and staining agents. All plates come with lids.

**ADDITIONAL INFORMATION**

<b>Weight</b>	3.83 kg
<b>Dimensions</b>	19.5 × 13.3 × 13.0 cm
<b>Unidades</b>	pore size 0.1µm (48 units), pore size 1.0 µm (48 units)

**THERE ARE NO REVIEWS YET.**

