

POLYESTER (PET) MEMBRANE TISSUE CULTURE PLATE INSERTS, 6 WELLS



42.26€ - 44.16€

Semitransparent polyester membrane with 6 wells and a pore size 0.1 μm , 0.4 μm , 1.0 μm , 3.0 μm , 8.0 μm or 12.0 μm . Sterilized. Packaging: 24 units (box with 4 packs of 6 units)

SKU: TCS010006 | **Category:** [Tissue Culture Plate Insert](#) |

VARIATIONS






Image	SKU	Price	Description	Unidades
	TCS010006	42.26€		pore size 0.1 µm (24 units)
	TCS011006	44.16€		pore size 0.4 µm (24 units)
	TCS012006	44.16€		pore size 3 µm (24 units)
	TCS013006	44.16€		pore size 8.0 µm (24 units)
	TCS015006	44.16€		pore size 1.0 µm (24 units)

Image	SKU	Price	Description	Unidades
	TCS100006	44.16€		pore size 12.0 µm (24 units)

GALLERY IMAGES



PRODUCT DESCRIPTION

Tissue Culture Plate Inserts feature a thin, microscopically polyester membrane that is tissue culture treated for optimal cell attachment and growth. Tissue Culture Plate Inserts provide excellent cell visibility under phase contrast microscopy and allow assessment of cell viability and monolayer formation. Tissue Culture Plate Inserts are sterile and is assembled with 6 well plates and a pore size 0.1 µm, 0.4 µm, 1.0 µm, 3,0µm, 8.0 µm or 12.0 µm. All plates come with lids.

ADDITIONAL INFORMATION

Weight	1.94 kg
Dimensions	19.5 × 13.3 × 13.0 cm
Unidades	pore size 0.1 µm (24 units), pore size 0.4 µm (24 units), pore size 1.0 µm (24 units), pore size 12.0 µm (24 units), pore size 3 µm (24 units), pore size 8.0 µm (24 units)

THERE ARE NO REVIEWS YET.

