

## ROLLER BOTTLE STANDARD, SURFACE TREATED, 1000 ML



**113.75€**

Roller bottle standard, surface treated and vent cap. Volume 1000 ml. Sterile. 24 units/box

**SKU:** TCB012001 | **Category:** [Roller Bottles](#) |

### PRODUCT DESCRIPTION

Roller Bottles are optimal applications for large yield of cell growth and multiplication, that can be used in both research and manufacturing, including the growing of mammalian cells for purposes of virus propagation and bio-product production.

Volume 1000 ml and wide cell growth area of 490 cm<sup>2</sup> for adhesions and huge growth volume for suspension.

Available surface treated.

Designed in large knurls on the cap for easy grip.

Vented cap style.

Bottle surface is smooth and free from striation to maximize usable area for growth.

The bottle has silk-printed graduation. Suitable for all the common instruments and automation.

1



