


## ROTATABLE CELL SCRAPER STERILE, 23 CM.



**52.30€**

Cell scrape with a rotatable cupboard, length 23 cm. 150 units/box individually packed.

**SKU:** CSC211023 | **Category:** [Cell Scraper](#) |

### PRODUCT DESCRIPTION

Rotatable cell scraper.

Changes of the blade angle of the cell scraper require slight pressure on the handle using the forfinger, thus pushing the handle downward towards the floor of the container.

Length of 23cm

Material: blades/PE, handle/ABS

Free rotating blade to twist to the desired direction

Total access to all corner  
Small raised knobs on the handle  
Individually wrapped  
Sterilized by gamma irradiation  
DNase/RNase-free & Non-pyrogenic

ADDITIONAL INFORMATION

Weight	1.80 kg
Dimensions	40.9 × 31.2 × 14.9 cm
Unidades	



