

## SEROLOGICAL PIPET, GRADUATED. INDIVIDUALLY PACKAGE



**16.91€ - 44.38€**

Serological Pipets Individually Package (paper/plastic). Material PS and Sterile. Volumes: 1, 2, 5, 10 and 25 ml. Standard or stretched mouth. Graduations: 1/100 ml, 1/50 ml, 1/10 ml, 2/10 ml, 5/10 ml. Packaging: 500 units (1 and 2 ml), 200 units (5 and 10 ml), 150 units (25 ml) and 100 units (50 ml)

**SKU:** GSP010010 | **Category:** [Serological Pipets](#) |

## VARIATIONS


Image	SKU	Price	Description	Unidades
	GSP010001	30.19€		1 ml (1/100 ml) Yellow (500 units)
	GSP010005	16.91€		5 ml (1/10 ml) Blue (200 units)
	GSP010010	21.74€		10 ml (1/10 ml) Orange (200 units)
	GSP010025	21.74€		25 ml (2/10 ml) Red (150 units)
	GSP010050	30.19€		50 ml (5/10 ml) Purple (100 units)

Image	SKU	Price	Description	Unidades
	GSP010102	44.38€		2 ml (1/100 ml) Green (500 units)
	GSP010002	36.23€		2 ml (1/50 ml) Green (500 units)
	GSP211010	19.32€		10 ml Stretched Mouth (1/10 ml) Orange (200 units)

GALLERY IMAGES



PRODUCT DESCRIPTION

Serological Pipets Individually Package (paper/plastic). Material PS, Sterile.  
Volumes: 1, 2, 5, 10 and 25 ml. Standard or stertched mouth  
Graduations: 1/100 ml, 1/50 ml, 1/10 ml, 2/10 ml, 5/10 ml.

Graduations are calibrated for accurate dispensing to within  $\pm 1\%$ .  
Different color ring for easy identification.  
Bidirectional graduations on the pipets provide added applicability.  
Negative graduation allows additional working volume.  
All pipets are supplied with a filter plug.  
Strict consistency tested.  
Non-pyrogenic.

ADDITIONAL INFORMATION

Weight	2.0 kg
Dimensions	39.8 × 31.2 × 11.0 cm
Unidades	1 ml (1/100 ml) Yellow (500 units), 10 ml (1/10 ml) Orange (200 units), 10 ml Stretched Mouth (1/10 ml) Orange (200 units), 2 ml (1/100 ml) Green (500 units), 2 ml (1/50 ml) Green (500 units), 25 ml (2/10 ml) Red (150 units), 5 ml (1/10 ml) Blue (200 units), 50 ml (5/10 ml) Purple (100 units)



