

SYRINGE-DRIVEN FILTERS MEMBRANE CA, SIZE 0,22MM. INDIVIDUALLY PACKED, STERILIZED



21.68€

Syringe-driven Filters membrane CA, size 0,22µm. Individually Packed, Sterilized. Packaging: 45 units

SKU: FCA206030 | **Category:** [Syringe Driven Filters](#) |

PRODUCT DESCRIPTION

Syringe driven filters are widely used in laboratory, which are easy to operate, safe and effective.

Membrane type CA
Pore size 0.22µm
Quality tested for filter efficacy and housing integrity
Low hold-up volume
Gamma radiation sterilized
Non-pyrogenic

Syringe driven filters are widely used in laboratory, which are easy to operate, safe and effective.

Quality tested for filter efficacy and housing integrity
Product Description Order Information.

1

THERE ARE NO REVIEWS YET.

