

SYRINGE-DRIVEN FILTERS, NYLON, 30 MM 0,22 MM, STERILE



21.29€

Syringe-driven Filters membrane Nylon, size 0,22µm. Housing diameter 30 mm. Sterilized.
Packaging: 50 units in bulk

SKU: FNY212030 | **Category:** [Syringe Driven Filters](#) |

PRODUCT DESCRIPTION

Syringe driven filters are widely used in laboratory, which are easy to operate, safe and effective.
Membrane type NYLON

30 mm and pore size 0.22µm
Quality tested for filter efficacy and housing integrity
Low hold-up volume
Sterilized
Non-pyrogenic

Syringe driven filters are widely used in laboratory, which are easy to operate, safe and effective.
Quality tested for filter efficacy and housing integrity
Product Description Order Information.

1

THERE ARE NO REVIEWS YET.

