

## SYRINGE-DRIVEN FILTERS, NYLON, 30 MM 0,22 MM, STERILE



**449.00€**

Syringe-driven Filters membrane Nylon, size 0,22µm. Housing diameter 30 mm. Sterilized.  
Packaging: 10 bags of 50 units in bulk

**SKU:** FNY212030 | **Category:** [Syringe Driven Filters](#) |

### PRODUCT DESCRIPTION

Syringe driven filters are widely used in laboratory, which are easy to operate, safe and effective.  
Membrane type NYLON

30 mm and pore size 0.22µm  
Quality tested for filter efficacy and housing integrity  
Low hold-up volume  
Sterilized  
Non-pyrogenic

Syringe driven filters are widely used in laboratory, which are easy to operate, safe and effective.  
Quality tested for filter efficacy and housing integrity  
Product Description Order Information.

1

**THERE ARE NO REVIEWS YET.**

