

SYRINGE-DRIVEN FILTERS, PES, 0,45 MM, STERILE



21.28€ - 24.09€

Syringe-driven Filters membrane PES, size 0,45µm. Housing diameter 13 or 30 mm. Bulk, Sterilized. Packaging available: 45 or 100 units

SKU: FPE404013 | **Category:** [Syringe Driven Filters](#) |

VARIATIONS

| Image | SKU | Price | Description | Unidades |
|---|-----------|--------|-------------|-------------------|
|  | FPE404013 | 24.09€ | | 13 mm (100 units) |
|  | FPE404030 | 21.28€ | | 30 mm (45 units) |

GALLERY IMAGES



PRODUCT DESCRIPTION

Syringe driven filters are widely used in laboratory, which are easy to operate, safe and effective.
Membrane type PES
13 or 30 mm and pore size 0.45µm
Quality tested for filter efficacy and housing integrity

Low hold-up volume
Sterilized
Non-pyrogenic

Syringe driven filters are widely used in laboratory, which are easy to operate, safe and effective.
Quality tested for filter efficacy and housing integrity
Product Description Order Information.

ADDITIONAL INFORMATION

| | |
|------------|-------------------------------------|
| Weight | 1.90 kg |
| Dimensions | 44.5 × 31.5 × 17.0 cm |
| Unidades | 13 mm (100 units), 30 mm (45 units) |

