

CENTRIFUGE TUBE, CONICAL BOTTOM AND EASYFLIP TOP CAP, 50 ML, STERILIZED



100.19€

Centrifuge tube conical bottom and EasyFlip top cap. Volume 50 ml. Sterilized, DNase/RNase free and Non-pyrogenic. Paper Rack

SKU: CFT212500 | **Category:** [Centrifuge Tubes](#) |

PRODUCT DESCRIPTION

Centrifuge tube conical bottom, EasyFlip top cap and sterilized. Volume 50 ml.
Made from high grade gamma resistance polypropylene
One hand easy flip to open the cap

Max.RCF: 9,400xg

Black printed graduations and white writing area can withstand chloroform

Autoclavable at 121°C and freezable to -80°C

Packaged in zip sealed bags or environmental friendly, recyclable paper racks

Leak proof

Gamma irradiation sterilized

Non-pyrogenic & DNase/RNase-free

ADDITIONAL INFORMATION

Weight	0.1 kg
Unidades	15 ml (Paper rack), 15 ml (re-sealable bag), 50 ml (re-sealable bag)

