

## MICRO CENTRIFUGE TUBE, CONICAL, 0,5 ML, STERILIZED







**5.44€ - 60.38€**

Micro Centrifuge tube conical bottom. Volume 0,5 ml. Sterilized, DNase/RNase free and Non-pyrogenic

**SKU:** CFT-001-005 | **Category:** [Micro Centrifuge Tubes](#) |

## VARIATIONS

Image	SKU	Price	Description	Unidades
	CFT-001-005	10.88€		0,5 ml (1000 units)
	CFT001015	5.44€		1,5 ml (500 units)
	CFT001020	6.53€		2 ml (500 units)
	CFT001050	60.38€		5 ml (1800 units)

**PRODUCT DESCRIPTION**

Micro Centrifuge tube conical bottom. Volume 0,5 ml. Sterilized.  
Engraved graduation ensure accuracy.  
Flat and frost caps surface together with smooth anf frosted body surface.  
DNase/RNase free and non-pyrogenic.

Autoclavable at 121°C and freezable to -80°C

ADDITIONAL INFORMATION

Weight	N/A
Dimensions	N/A
Unidades	0,5 ml (1000 units), 1,5 ml (500 units), 2 ml (500 units), 5 ml (1800 units)



