

ULTRA FREEZER -86°C LABCARE RANGE, 700 LITERS



7,967.00€

Ultra-freezer -86°C with VIP panels and double refrigeration circuit technology with autonomous cascade compressor that allows it to work at -70°C in case of failure of the other circuit. Without oil separator. Capacity 700 Liters.

SKU: ULF70086 | **Category:** [Equipments](#) |

PRODUCT DESCRIPTION

Ultra-freezer -86°C with VIP panels and double refrigeration circuit technology with autonomous

cascade compressor that allows it to work at -70°C in case of failure of the other circuit. Without oil separator.

Digital Controller with Touch Screen that integrates a temperature logger with graph registration and display, data export via USB, access control by password, and setting of visual and audible alarms.

Door with double perimeter magnetic seal and heated.

Main Features:

- 2 Bulwarks
- Front document holder
- 4 front doors and 3 height adjustable shelves.
 - Noise 56dBA
 - Ergonomic lock handle.
 - Vacuum break point
- Temperature range from -40 ° C to -86 ° C with Datalogger, graphic representation and USB to facilitate data collection.
- Maximum protection thanks to visual and audible alarms due to power failure, drop or rise in temperatures, probe failure, open door, filter saturation and discharged battery.
 - Blind door with lock and 4 interior doors to avoid cold losses.
- Audited and certified by an international organization AENOR in ISO 9001: 2008 (Quality); ISO 14001: 2004 (Environment); OSHAS 18001: 2007 (Safety and Health for workers).
 - Storage capacity cryoboxes of 2'' (5cm): 576 units

Technical specifications:

- SPECTRUM LABCARE
 - MODEL ULF70086
 - RANGE Tª 0-40°C / -86°C
- NOMINAL CAPACITY (Liters) 700
- GROSS CAPACITY (Liters) 778
- NET CAPACITY (Liters) 751
- EXTERIOR MEASURES (mm) 1224 x 915 x 1986
- INTERIOR MEASURES (mm) 920 x 650 x 1300
- DOOR CLEAR (mm) 920 x 1297
 - N ° SHELVES 3
 - N ° OF LEVELS 4

1

THERE ARE NO REVIEWS YET.

