

VACUUM-DRIVEN BOTTLES, FCA 0,22MM, 250 ML.



38.04€

Vacuum-driven Bottles, membrane FCA, pore size 0,22µm and volume 250 ml. Packaging: 12 units

SKU: FCA206250 | **Category:** [Vacuum Filtration](#) |

GALLERY IMAGES



PRODUCT DESCRIPTION

Vacuum-driven Bottles, membrane FCA, pore size 0,22µm and volume 250 ml. Membrane pore size 0.22µm

Light weight and heavy wall construction

Large knurls on the reservoir bottle cap for easy screw

Designed wide and easy access bottle mouth for efficiently and stably pour out

Engraved graduation ensure veracity

Ergonomically designed sidewalls and collar can simplify tightening/loosening and adjustments

Designed hose connector can fit multiplicate hose diameters

Non-pyrogenic

1

