

ZEROTIP PIPPETE MICRO TIP WITHOUT FILTER 10 ML, NON-STERILE



80.60€

ZEROTIP Pipette micro tip low retention without filter 10 µl and Non-sterile. Natural. Packaging:
2 boxes with 10 racks x 96 units

SKU: PMT250010 | **Category:** [ZEROTIP Pipette micro tips](#) |

PRODUCT DESCRIPTION

ZEROTIP Pipette micro tips without filter 10 µl and Non-sterile. Natural. Rack box.

JET ZEROTIP® are manufactured from super clear high quality Polypropylene. The surfaces of the tips

are produced through a special process. This process makes the tip inner surface become super hydrophobic, thus significantly reduces sample loss and provides substantially higher reproducibility when working with critical media.

Free of Nucleic Acids, Pyrogens/Endotoxins, PCR Inhibitors and Trace Metals.

Ideal for biological samples that contain detergents such as SDS, Tween Triton® -100 ect.

Tips with PP filter are also available.

Preferred accessories for most brand micropipettor.

Non-Pyrogenic & DNase/RNase-free.

1

THERE ARE NO REVIEWS YET.

