

## ZEROTIP PIPPETE MICRO TIP WITHOUT FILTER, STERILE





**80.60€**

ZEROTIP Pipette micro tip low retention without filter 10 µl and sterile. Natural. Packaging: 2 boxes with 10 racks x 96 units

**SKU:** PMT251010 | **Category:** [ZEROTIP Pipette micro tips](#) |

## VARIATIONS

Image	SKU	Price	Description	Unidades
	PMT251010	80.60€		10 µl (20 racks x 96units)
	PMT251200	80.60€		200 µl (20 racks x 96units)

GALLERY IMAGES



PRODUCT DESCRIPTION

ZEROTIP Pippete micro tips without filter 10 µl and sterile. Natural. Rack box.  
JET ZEROTIP® are manufactured from super clear high quality Polypropylene. The surfaces of the tips are produced through a special process. This process makes the tip inner surface become super

hydrophobic, thus significantly reduces sample loss and provides substantially higher reproducibility when working with critical media.

Free of Nucleic Acids, Pyrogens/Endotoxins, PCR Inhibitors and Trace Metals.

Ideal for biological samples that contain detergents such as SDS, Tween Triton® -100 ect.

Tips with PP filter are also available.

Preferred accessories for most brand micropipettor.

Non-Pyrogenic & DNase/RNase-free.

**ADDITIONAL INFORMATION**

<b>Weight</b>	2.30 kg
<b>Dimensions</b>	0.30 × 0.24 × 0.26 cm
<b>Unidades</b>	10 µl (20 racks x 96units), 1000 ul (2 boxes of 10 Racks x 96 units), 1000 ul long (2 boxes of 10 Racks x 96 units), 200 µl (20 racks x 96units), 300 ul (2 boxes of 10 Racks x 96 units)



